Phone: +91-8860605265 Email: info@astrascientific.com

**Product Name:** 

Astra Scientific Free and Forced Vortex

**Product Code:** 

Fluid Mechanics Lab0016



## **Description:**

Astra Scientific Free and Forced Vortex

## **Technical Specification:**

## **DESCRIPTION:**

The Free and Forced Vortex Apparatus for the study of the shape of 'free and forced vortices' consists of a cylindrical, transparent vessel, having two pairs of diametrically opposed inlet tubes of 9.0 mm and 12.5 mm diameter. The 12.5 mm diameter inlet tubes which are angled at 15° to he diameter, so that a swirling motion is imparted to the liquid entering the vessel, are used as entry tubes for the free vortex experiment.

A smooth outlet is centrally positioned in the base of the vessel and a set of push-in orifices of 24, 16, 12 and 8 mm diameter is supplied to reduce the outlet diameter to a suitable value. A gauge determines the profile of the vortex formed at the top of the vessel, housed on a diametrically mounted bridge piece, which measures the diameter of the vortex at various depths. This gives the co-ordinate points required to plot the vortex profile.

Phone: +91-8860605265

Email: info@astrascientific.com

The forced vortex is created in the vessel described above by using as the inlet the 9 mm bore tubes which are angled at 60° to the diameter. The input water from these tubes impinges on a simple two-blade paddle, which acts as a stirrer/flow straightener. The water 'leaves' the vessel via the 12.5 mm diameter angled tubes, which are used as the 'entry' tubes for the free vortex experiment.

#### **FEATURES:**

- Experiment to plot the shape of a free vortex by measurement of the surface profile co-ordinates, and thus verify that v r = constant where v is the speed and r is the radius of the vortex
- Experiment to plot the surface profiles of various forced vortices formed under different speed conditions
- Verification of the formula h=(w² r²)/2g for forced vortices where h
  is the height of the surface of the water above the datum point; w is
  the vortex angular velocity and r is the vortex radius.

#### SPECIFICATIONS:

- Cylindrical, transparent vessel 250 mm diameter &180 mm depth,
- Inlet tubes of 9.0 mm and 12.5 mm diameter.
- Set of push-in orifices of 24, 16, 12 and 8 mm diameter



# **Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265