Phone: +91-8860605265 Email: info@astrascientific.com

Product Name :Monsanto Hardness Tester

Product Code:
PHARMACEUTICAL SUKAS03LAB012A



Description:

Monsanto Hardness Tester

Technical Specification:

This range is appreciated for applications in numerous industries and sectors for various purposes. Hardness testing is a laboratory technique used by the pharmaceutical industry to test the breaking point and structural integrity of a tablet "under conditions of storage, transportation, and handling before usage" The breaking point of a tablet is based on its shape. Quality analysts at our end make sure to fabricate this range by following quality standards and guidelines in sync with the industry set norms.

Phone: +91-8860605265

Email: info@astrascientific.com

_					
-	ea	tı	ır	0	2 .

Adjust the Scale (by sliding), so that the zero on the Scale coincides with the pointer.

Turn the screw knob slowly till the tablet breaks.

Hold the Tablet between the jaw and nozzle in edgewise position.

The pressure indicated on the dial is in Kilogram per sq. cm.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265