

Product Name :
Sterility Test Apparatus

Product Code :
PHARMACEUTICAL SUK-
AS03LAB097A



Description :

Sterility Test Apparatus

Technical Specification :

The liquid sample is drawn through this sterile membrane filter to concentrate any contaminating micro-organisms on the filter surface. The sterility tests in these systems are done by the bacteria-retentive membranes filter method. A 47 mm membrane filter is installed in a suitable holder. After filtration is completed, sterile flushing fluid is passed into the funnels, rinsing all sample residues through the filters so that it cannot possibly inhibit growth of the retained micro-organisms. After this, the filter is sectioned and aseptically transferred to media for culturing. If growth is not visible after the incubation period, the sample passes the test.

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific