

**Product Name :**  
Standard Vicat Apparatus

**Product Code :**  
ASTCEMENT-0020



**Description :**

Standard Vicat Apparatus

**Technical Specification :**

**General Description:** The

Vicat frame consists essentially of a metal stand with a sliding rod. An adjustable indicator moves over a graduated scale. The needle or plunger is attached to the bottom end of the rod to make up the test weight of 300 g.

The frame is supplied with accessories, depending on the requirement. It is also proposed with basic EN accessories or ASTM/AASHTO accessories. Weight approx.: 4 kg

Vicat test set, EN method, including mould, 1.13 mm dia. needle, consistency plunger 10 mm dia., supporting plate, glass thermometer and final needle.

Vicat test set, ASTM/AASHTO method, including mould, 1 mm dia. needle, consistency plunger 10 mm dia., supporting plate and glass thermometer.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265