


Product Name :
HYDROSTATIC PRESSURE, Tilting Submerged Quadrant

Product Code :
H6



Description :

Quadrant inner radius 100 mm, outer radius 200 mm, and width 75 mm. Tilting up to 30 degrees in either direction. Balance arm with counter weight, and weights. Transparent tank 8 lt

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

