

Product Name :
Astra Scientific Digital Turbidity Meter

Product Code :
ASTENVIOR-0007



Description :

Astra Scientific Digital Turbidity Meter

Technical Specification :

DIGITAL TURBIDITY METER

DISPLAY 3.5 Digital LED.

POWER 230v $\hat{\pm}$ 10% 50 Hz.

RANGE 0 to 200 NTU.

RESOLUTION 1 NTU.

ACCURACY $\hat{\pm}$ 2% FS $\hat{\pm}$ 1 Digit.

SAMPLE SYSTEM 30mm Clear Glass Test Tubes.

LIGHT SOURCE 6V-1Amp. Tungsten Lamp.

DETECTOR Photocell / Photodiode.

DIGITAL TURBIDITY METER

TURBIDITY RANGE 0 -1 NTU, 0-10 NTU,

0-100 NTU, 0-1000 NTU

ACCURACY UP To $\pm 3\%$,

DETECTOR PHOTO DIODE

BULB 6 V , 6 W

DISPLAY 3 $\frac{1}{2}$ DIGIT PANEL METER

POWER 220 volts, 50 Hz AC

ACCESSORIES 1) TEST TUBE 4 NOS

2) CELL RISER

3) LIGHT SHIELD

4) INSTRUCTION M



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265