

Product Name :
Astra Scientific Digital Dissolved Oxygen

Product Code :
ASTENVIOR-0008



Description :

Astra Scientific Digital Dissolved Oxygen

Technical Specification :

DIGITAL DISSOLVED OXYGEN ANALYZER CUM TEMP. METER

RANGE 0 to 20.0ppm

RESOLUTION 0.1 ppm

ACCURACY ± 0.2 ppm

TEMPERATURE COMP. 0to 50 deg. C. Manual.

DISPLAY 3.5 Digit LED

STIRRER Provided.

POWER 230V $\pm 10\%$ 50 Hz

SENSOR DO Amperometric Gold / Silver Temp. RTD (PT-100)

TIME RESPONSE 10 Sec

DIGITAL DISSOLVED OXYGEN METER

DISPLAY 3.5 Digit LED

RANGE 0 to 20.0ppm.

RESOLUTION 0.1 ppm.

ACCURACY $\hat{\pm}$ 0.2ppm

TEMPERATURE COMP. 0 to 50 deg. C.

SENSOR DO Amperometric Gold / Sliver

TIME RESPONSE 10 Sec.

POWER 230V $\hat{\pm}$ 10% 50 Hz AC



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265