

Product Name :
Astra Scientific Cantilever Beam Apparatus

Product Code :
Applied-Lab-0004



Description :

Astra Scientific Cantilever Beam Apparatus

Technical Specification :

Material	Stainless Steel
Weight set	100 gm to 1 kg
Horizontal arm	75 cm
Vertical arm	20 cm
Display	Analog
Frequency	50 Hz

Consists of a Stainless steel beam of L type shape with horizontal graduated horizontal arm 75cm and vertical arm of 20cm, with a spring balance of 10 kg capacity to make the horizontal leg horizontal under any load. Complete with one sliding slotted set weight set 100gm to 1 kg.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265