

Product Name : Astra Scientific Section Model of Two Stroke And Four-Stroke Petrol Engine	Product Code : AutomobileLab0022
---	--



Description :

Astra Scientific Section Model of Two Stroke And Four-Stroke Petrol Engine

Technical Specification :

The section model represents a real 2-stroke and 4-stroke petrol engine, cut in every detail, carburetor, ignition, etc. The suction-exhaust-transfer channels are especially highlighted so as to make it easy to learn the cycle.

This cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc. Many parts have been chromium-plated and galvanized for a longer life.

FEATURES:

- Operated: Manually (Through a Crank Handle)
- Equipment Layout: On Base
- Fine finish

- Simple operations
- Functional accuracy

SPECIFICATIONS:

- TWO STROKE CUTSECTION ENGINE
 - o Piston displacement 48 cc
 - o Air cooling
 - o Plug point and magnet flywheel ignition
 - o Box carburettor
- FOUR STROKE CUTSECTION ENGINE
 - o Displacement: 160cc
 - o Camshaft in the crankcase
 - o Overhead valves
 - o Rpm regulator
 - o Oil pump
 - o Carburetor
 - o Air filter
 - o Silencer
 - o Tank



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265