

**Product Name :**  
Astra Scientific CNC Plasma Cutting Machine

**Product Code :**  
CADLab0002



**Description :**

Astra Scientific CNC Plasma Cutting Machine

**Technical Specification :**

The advanced mechanism designs ensure the high performance of the dynamic response. Reduce the finishing work- as the cutting error is very small, so can reduce the work piece reworking. Automatic height tracking system- Automatic level control systems with fast response characteristics have guaranteed the quality under the condition of high-speed...

**FEATURE**

- High Stability- the advanced mechanism designs ensure the high performance of the dynamic response.
- Reduce the finishing work- as the cutting error is very small, so can reduce the work piece reworking.
- Automatic height tracking system- Automatic level control systems with fast response characteristics have guaranteed the quality under the condition of high-speed cutting.
- Easy Installation and maintenance- Reasonable modular design is convenient for assembly and maintenance.

SPECIFICATION:

- Max Cutting Length: 0-500 mm, 500-1000 mm
- Automation Grade: Automatic
- Machine Type: CNC Cutting Machine
- Finishing: Powder Coated
- Max Cutting Speed: 500-1000 mm/min
- Max Cutting Thickness: 15 mm
- Max Cutting Width: 500-1000 mm
- Accuracy:  $\pm 0.02$ mm
- Repeatability precision:  $\pm 0.02$ mm



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265