


Product Name :
Astra Scientific CNC Plasma Cutting Machine

Product Code :
CADLab0005



Description :

Astra Scientific CNC Plasma Cutting Machine

Technical Specification :

The advanced mechanism designs ensure the high performance of the dynamic response. Reduce the finishing work- as the cutting error is very small, so can reduce the work piece reworking. Automatic height tracking system- Automatic level control systems with fast response characteristics have guaranteed the quality under the condition of high-speed...

FEATURE:

- High Stability- the advanced mechanism designs ensure the high performance of the dynamic response.
- Reduce the finishing work- as the cutting error is very small, so can reduce the work piece reworking.
- Automatic height tracking system- Automatic level control systems with fast response characteristics have guaranteed the quality under the condition of high-speed cutting.
- Easy Installation and maintenance- Reasonable modular design is convenient for assembly and maintenance.

SPECIFICATION:

- Max Cutting Length: 0-500 mm, 500-1000 mm
- Automation Grade: Automatic
- Machine Type: CNC Cutting Machine
- Finishing: Powder Coated
- Max Cutting Speed: 500-1000 mm/min
- Max Cutting Thickness: 15 mm
- Max Cutting Width: 500-1000 mm
- Accuracy: ± 0.02 mm
- Repeatability precision: ± 0.02 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265