

Product Name :
Astra Scientific Super Scanner

Product Code :
Civil0002



Description :

Astra Scientific Super Scanner

Technical Specification :

Handheld Metal Detector

1. For fast weapon screening and theft prevention
2. Highly sensitive scanning circuit-fast scan rate
3. Instant response-clear visual and audio alarms
4. No false alarms-digital go, no-go threshold circuit
5. Extra large coil for wide coverage
6. High sensitivity for non-intrusive body searching
7. Ergonomically designed-reduces operator fatigue
8. Electronic switching with tactile feedback
9. Three state, daylight-visible, Blue LED alarm and status indicator

10. Low field strength-no effect on pacemakers and magnetic media

11. Detects all metals, alloys and ferrites

12. Low current drain-over 300 hrs of use on 9V Duracell battery



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265