

Product Name :
Astra Scientific Digital Multimeter

Product Code :
Digital0001



Description :

Astra Scientific Digital Multimeter

Technical Specification :

Digital Multimeter

We are Leading Supplier of Digital Multimeter in India. We also bulk supply Digital Multimeter for laboratory tensors.

Features

- 3 1/2 Digit Manual Ranging
- Temperature Measurement
- Capacitance Measurement
- Data hold and auto power off function

Technical Specifications :

DC Voltage :

Range : 200mV ~ 600V

Accuracy : (0.5% + 1digit)

AC Voltage :

Range : 2V ~ 600V

Accuracy : (0.8% + 3digits)

DC Current :

Range : 2mA ~ 10A

Accuracy : ;(0.8% + 1digit)

AC Current :

Range : 20mA ~ 10A

Accuracy : $\pm(1.2\% + 5\text{digits})$
Resistance :
Range : 200 ~ 20M Ω
Accuracy : (1.0% + 2digit)
Diode test open circuit : 2.8 volt
Capacitance :
Range : 2nF ~ 200 μ F
Accuracy : (4.0% + 3digits)
Temperature :
Range : -40 $^{\circ}$ C ~ 1000 $^{\circ}$ C
Accuracy : $\pm(1\% + 7\text{digits})$
Power : DC 9 V battery



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265