

**Product Name :**  
Resistance Box

**Product Code :**  
PHYSICSSUK-AS02LAB004



**Description :**

Resistance Box

**Technical Specification :**

Resistance Box that is designed to perform multiple operations. Owing to its features, this precision instruments cover the requirements of both field and laboratory work. This box is designed using optimum quality components & cutting-edge techniques in accordance with industry standard. This resistance box is designed with multiple plugs for measuring the electrical values like voltage, current and resistance in a hassle-free and accurate manner.

**Range:** 0.1 to 1, 1 to 10, 10 to 100, 100 to 1000, 1000 to 10000 OHMS.

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**