

**Product Name :**  
Vernier Calliper

**Product Code :**  
PHYSICSSUK-AS02LAB005



**Description :**

Vernier Calliper

**Technical Specification :**

Designed using the best quality raw material, this caliper is used in Physics Lab and Science labs. It is a precision instrument that can be used to measure internal and external distances extremely accurately.

**Features:-**

- Length: 12.5cm
- Stainless Steel

- Resistant to abrasion
- Accurate measurement
- Robust construction
- Smooth finish.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265