

Product Name :
Hand Grinding Mill

Product Code :
PHARMACEUTICAL SUK-
AS03LAB006A



Description :

Hand Grinding Mill

Technical Specification :

Makes Fresh and Healthy milled grain products in minutes. Complete on metallic stand having adjustable arrangement for making coarse & fine powder. Body Made of very strong cast iron. Secures to your table or counter top. Highly recommended for preparation of samples in pharmacy labs.

Features:-

Easy turn handle is made of hardwood.

Easy to clean and assemble

Tall Large Head Cup

Perfect grinds corn/coffee/grains.

Adjustment knob on side to adjust from coarse to fine



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265