

**Product Name :**  
Tar Viscometer Automatic

**Product Code :**  
FUEL SUK-AS03LAB005



**Description :**

Tar Viscometer Automatic

**Technical Specification :**

Tar Viscometer Automatic: The apparatus consist of bath with cup of 10 mm or 4 mm orifice and sleeve stirrer with ball lifting clip and ball .The bath is fitted with an immersion heater to take the water to the required temperature and a drain valve. The temperature is controlled by energy regulator



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265