

Product Name :
Colour Matching Cabinet

Product Code :
TEXTILE SUK-AS03LAB004



Description :

Colour Matching Cabinet

Technical Specification :

For visual colour matching of yarn, fabric and garments under different imported light sources like D-65, UV, TL-84, CWF, TL-83, U-30, U-35, Inca A, with Electronic ballast & hour meter metal body without door available standards : ASTM D 1792, ISO 3664, BS 950



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265