

Product Name :
Fully Automatic Twist Tester Micro Processor Based

Product Code :
YARN SUK-AS03LAB001



Description :

Fully Automatic Twist Tester Micro Processor Based

Technical Specification :

Available for both Single as well as double yarn.

As per standards:- ASTM D-1422/1423,

ISO-2061, BS-2085



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265