

Product Name :
Ross Flex Tester

Product Code :
FOOTWEAR SUK-AS03LAB004



Description :

Ross Flex Tester

Technical Specification :

The ability of shoes or chappals to withstand repeated bending experience during walking without developing a crack within their estimated life span is a very important factor for determining the quality of its sole. Also, the ability of the sole to withstand propagation of any slit formed due to a sharp edge or pointed object during walking is equally important in assessing its quality.

Technical specification

Bend Radius

20 Age
50 Hz
150 Specimen Diameter
3.5 Width
1000 Cycles Per Minutes
90 Degree Specimen Angle



Astra Scientific

Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265