

Product Name :
Bennewart Flex Tester

Product Code :
FOOTWEAR SUK-AS03LAB007



Description :

Bennewart Flex Tester

Technical Specification :

Bennewart Flex Tester.

These equipments are manufactured by our highly skilled manpower using sophisticated technology. Our quality analyzer stringently checks these flex testing equipments on quality parameters before delivering to our clients. **Bennewart Flex Tester** is broadly used to measure the resistance of a component to cut growth during repeated flexing. Clients can buy these equipments from us t cost effective prices.

Further details:

It can also be used to assess the effect of surface patterns on crack initiation and growth. This test is especially applicable to the outsoles of footwear, but may also be used with certain other flexible components.

Technical specification

Bending Mandrel Diameter
Display Type
Material
Grade Automatic
Voltage
Frequency
Brand Enterprises
Operation Frequency

File Name
Digit



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265