

Product Name :
Crimp Tester

Product Code :
YARN SUK-AS03LAB006



Description :

Crimp Tester

Technical Specification :

With great support from our talented employees, we have been engaged in the business of providing Crimp Tester.

It is utilized for testing presence of small folds, curls, or undulations in fabric yarn strands. This tester is used by gripping a piece of specimen to clamps present on both sides of machine. Then tension is applied at preset rate with real time display of force and elongation value. Crimp Tester is greatly admired for having easy to use nature and long life.

Features:

- Equipped with touch screen control panel
- Auto stop and alarms are provided after predetermined drawing of tension value
- Compact and strong structure



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265