

Product Name :
Zoo Plankton Net

Product Code :
LABUK-AP12LAB0005



Description :

Zoo Plankton Net

Technical Specification :

- The upper part of the net is made of strong Beaver Nylon.
- Standard mounting with 3pcs of polyester line 03 mm. with bridle, length 60 cm.
- The conical net bag with an upper diameter of 30 cm and a length of 100 cm is made of Monodur Nyltal (Nylon), mesh size 20 to 200 μ ,m

- The sample can be removed through a valve at the bottom of the sample cylinder.
- Sample cylinder made of PVC, 5 windows covered by net, mesh size according to the plankton net, minimum mesh size of 20 μm .
- Content 80 ml.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265