

**Product Name :**  
Soil Ph Meter

**Product Code :**  
LABUK-AP12LAB0006



**Description :**

Soil Ph Meter

**Technical Specification :**

**Features:**

Temperature range of -10 to 122°F (-15 to 50°C) and is housed in a plast.

Range: pH: -2.0 to 16.0

pH: Temperature: -5 to 105oc

Resolution: pH: 0.01

Temperature : 0.1oc

Accuracy: pH @ 20C:  $\pm 0.02$

Temperature: 0.5 °C

Calibration: automatic one or two point

Calibration with two sets of memorized buffers (Standard 4.01, 7.01, 10.01).

Temperature Compensation: Automatic, - 5 to 105°C

Battery Type / Life: 1 .5V AAA (3) / 1200 hours of continuous use

Auto-off after 8 minutes of non-use

Environment: 0 to 50°C; RH max 100%

Dimensions: 152 x 58 x 30 mm

Weight: : 205 g

IC molded case and reservoir.



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265