


**Product Name :**  
Stopwatch Digital

**Product Code :**  
LABUK-AP12LAB0015



**Description :**

Stopwatch Digital

**Technical Specification :**

**Features:**

- 24-hour clock operation
- Stopwatch/Chronograph mode with 1/100sec resolution

- Programmable alarm with an hourly chime
- Calendar displays day, month, and date.
- Water resistant.
- Material: plastic.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265