

Product Name :
Laboratory Water Taps

Product Code :
LABUK-AP12LAB0019



Description :

Laboratory Water Taps

Technical Specification :

Features:

- Type: three way bench laboratory faucet.
- Material: metal thick solid brass.

- Coating: gloss epoxy powder, resistant to most chemicals, UV fading and Heat.
- Ceramic Core of Valves: 500000 life cycles, water static pressure Max. 25bar.
- Knob: high grade PP, humanistic design, convenient and shortcut.
- Nozzle: detachable nozzle, water flow control outlet attachable.
- Warranty for 2 years.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265