

**Product Name :**  
Bench Top Conductivity Meter

**Product Code :**  
LABUK-AP12LAB0029



**Description :**

Bench Top Conductivity Meter

**Technical Specification :**

**Features:**

- Conductivity range: 0 $\mu$ S/cm to 100  $\mu$ S/cm
- Accuracy:  $\pm 0.5$  % of value
- Temperature: 5.0 to 105.0oc  $\pm 0. 1$ oC (23 ... 22) °F)
- Power: 230V, 50Hz
- Display: digital LED
- Cell constant: 0.25 to 2.5 cm- 1. fix 0.01cm-1
- Calibration: 1-point

**Accessories:**

- Dust cover,
- Conductivity cell.
- Two conductivity standards of 12.88 mS/cm and 1413  $\mu$ S/cm

- Warranty: two year.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265