

**Product Name :**  
Pycnometer

**Product Code :**  
AST-GEO-0001



**Description :**

Pycnometer

**Technical Specification :**

Specific Gravity is the ratio of the weight in air of a given volume of a material at a standard temperature to the weight in air of an equal volume of distilled water at the same stated temperature. The Pycnometer Bottle can be used to test a wide range of materials from clay to sand and gravel smaller than 10mm.

It consists of 1kg Glass Jar with Brass Cone, Locking Ring and Rubber Seal.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265