

Product Name :
Blood Donor Chair Fixed Seat

Product Code :
HOSPITALSUK-AS120003



Description :

Blood Donor Chair Fixed Seat

Technical Specification :

Made of mild steel tubing

Padded seat, backrest and armrest

Fixed seat height

Adjustable armrest in any position for blood drawing

Maximum weight load 150 kgs

Without castors



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265