

Product Name :
Saline with I.V. Stand Double Hook

Product Code :
HOSPITALSUK-AS270001



Description :

Saline with I.V. Stand Double Hook

Technical Specification :

Height adjustable from 137 to 244cm

Strong tubular stand mounted on rectangular pipe with five castors of 5cm dia.

Double hook adjustable with a bakelite knob

Base made of PU, with stainless steel inner rod and outer pipe of metal



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265