


Product Name :
Saline with I.V. Stand Double Hook Rectangular Pipes

Product Code :
HOSPITALSUK-AS270002



Description :

Saline with I.V. Stand Double Hook Rectangular Pipes

Technical Specification :

Height adjustable from 137 to 244cm

Strong tubular stand mounted on rectangular pipes with five castors of 50cm dia.

Double hook adjustable with a bakelite knob

Base, inner and outer rod of complete stainless steel



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265