

Product Name :
X-Ray View Box

Product Code :
HOSPITALSUK-AS330001



Description :

X-Ray View Box

Technical Specification :

Used to view MRI, CT Scan and X Ray images

Energy efficient

Film active switch

Cost effective

Hospital X-Ray View Box - OEM Hospital Medical Hot Sale Super Thin X Ray Film Viewer Film Viewing Illuminator View Box with LCD Screen Digital Control Negatoscope FOB Price. We are a well-known manufacturer and supplier of X-Ray View Box which are required in hospitals, medical institutes and nursing homes. High-quality X-Ray View Box (LED) manufacturer and supplier from India. We offer a comprehensive range of X-Ray View Box at best prices. LED X Ray View Box, X Ray View Box manufacturers and OEM manufacturers India. Contact details of XRay View Box companies manufacturing and supplying company India China UK USA Dubai UAE



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265