

**Product Name :**  
Astra Gradient PCR

**Product Code :**  
ASTRALAB-SUPPLYA1011



**Description :**

Astra Gradient PCR

**Technical Specification :**

Amplification analyzer has strong comprehensive performance, fast amplification speed and long service life. It can also guarantee the reliability of the analysis results and the system when it is tested in more complex and harsh environment. The number of cycles is more than 200,000. This series of products have stable performance, simple operation, practicality and economy.

The most advanced peltier-based semiconductor technology.

High-performance universal power supply.

Large 5.7-inch high-definition LCD display.

Bilingual curve graphical interface.

Power-down data protection.

Metal shell, solid, practical, beautiful and generous.

Stepless adjustable Lid can adapt to different heights in vitro.

Lid can be positioned at any angle.

High sealing reaction zone to ensure stable and reliable experiment.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265