

**Product Name :**  
Astra Refrigerated Centrifuge

**Product Code :**  
ASTRALAB-SUPPLYA1014



**Description :**

Astra Refrigerated Centrifuge

**Technical Specification :**

A centrifuge is a machine that uses centrifugal force to accelerate the separation of different materials that need to be separated. Refrigerated centrifuges are divided into low-speed, high-speed refrigerated centrifuges, and ultra-speed analysis, preparation dual-use refrigerated centrifuges and other models.

With frequency conversion motor drive and microcomputer control

High-definition LCD and digital dual-screen display

With electronic safety door lock, imbalance protection

With automatic recognition of the rotor to avoid speeding

It has a unique elastic rotating head fixing method, which is convenient and quick to replace the rotating head

With 12 program storage, 40-level lifting speed is adjustable

With imported compressor, non-fluorine refrigerant, refrigeration/heating dual-loop control, precise temperature control

The RCF function can be directly set and displayed without RPM/RCF conversion

With imported brand-name frequency conversion controller, accurate speed regulation, durable

The operating parameters of the centrifuge can be modified during operation

The refrigerated centrifuge is suitable for high-speed angular rotors and low-speed horizontal rotors

All-steel body, three-layer protective steel sleeve to ensure safety



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265