

Product Name :
Astra Ultrasonic Homogenizer

Product Code :
ASTRALAB-SUPPLYA1016



Description :

Astra Ultrasonic Homogenizer

Technical Specification :

Ultrasonic homogenizers are ideal for protein and enzyme extraction.

The ultrasonic homogenizer JY92-IIN is an efficient and low-cost laboratory homogenizer.

With 650 W power, 20-25KHz converter, 6mm cone probe, it can crush and homogenize 0.5-500ml samples.

The instrument adopts new software and is controlled by a central computer

Adopt 7-inch high-definition TFT touch screen display

Ultrasonic power step, which can be adjusted continuously by 1%.

It has pulse and continuous operation functions, and has accurate test functions

Ultrasonic test time, gap time and total time can be accurate to 0.1 second.

Integrated anechoic chamber to reduce cavitation noise during processing

All D models with temperature indicator and controller

High quality reasonable price, time setting



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265