


**Product Name :**  
Astra PCR Plate

**Product Code :**  
ASTRALAB-SUPPLYA1019



**Description :**

Astra PCR Plate

**Technical Specification :**

This microplate qPCR plate is a high-quality type, thin wall can achieve the best heat transfer, for PCR research.

Widely used in genetics, biochemistry, immunity, medicine and other fields, in addition to basic research such as gene isolation, cloning and nucleic acid sequence analysis, it is used for disease diagnosis or any place with DNA and RNA, which is a disposable consumable in the laboratory. Product.

With high transparent PP material, thin and uniform wall to achieve the best heat transfer

Can have a clear number, easy to place

Provide silicone cover



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265