


Product Name :
Astra Mini Diaphragm Pump

Product Code :
ASTRALAB-SUPPLYA7009



Description :

Astra Mini Diaphragm Pump

Technical Specification :

The small diaphragm pump is a two-stage pump with medium gas. All parts in contact with gas are made of PTFE material. It has high corrosion resistance and a wide range of applications. It is suitable for chemical, pharmaceutical, petrochemical and other industries. Treatment of corrosive gases, such as suction filtration, vacuum distillation, rotary evaporation, vacuum concentration, centrifugal concentration, solid phase extraction, etc.

Vacuum is stable and easy to control

Fully resistant to chemical corrosion

No water consumption, strong environmental performance

Small and light

Easy to use and clean

Lasting and stable work

Easy to maintain and repair



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265