

Product Name :
Astra Industrial Peristaltic Pump

Product Code :
ASTRALAB-SUPPLYA7011



Description :

Astra Industrial Peristaltic Pump

Technical Specification :

The industrial peristaltic pump adopts a transparent protective cover, which can easily observe the operation of the pump head and hose. The flow can be controlled between 1.3~35L/Min. The speed is between 30rpm and 350rpm.

No pollution: the fluid only touches the pump tube, not the pump body;

High precision: good repeatability and high stability;

Good air tightness: it has good self-priming ability, can be left unused, and can prevent backflow;

Easy to maintain: no valves and seals;



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265