

Product Name :
Medical Thermostatic Incubator

Product Code :
ASTRALAB-SUPPLYA10001



Description :

Medical Thermostatic Incubator

Technical Specification :

Thermostatic Incubator is used in scientific research and industrial production departments such as medical and health, pharmaceutical industry, biochemistry and agricultural science for bacterial cultivation, fermentation and constant temperature testing.

Internal volume: 50L~600L;

Temperature control range: RT+5?65?, temperature fluctuation: $\pm 0.5?$; temperature uniformity: $\pm 1.5?$ (test point is 37?);

Controller: Dedicated temperature controller, fast induction, small system error

Data printing: Optional printer or RS485 interface; can connect to printer or computer;

Double-door structure: There is a built-in glass door for easy observation of culture changes;

Box material: The inner tank is made of brushed stainless steel, the outer shell is made of steel plate sprayed plastic, and the shelf can be freely loaded and unloaded;

Install power supply: AC 220V $\pm 10%$ 50HZ;

Working environment: +5 ~ 35?;



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265