

**Product Name :**  
Medical High Pressure Vessel

**Product Code :**  
ASTRALAB-SUPPLYA12005



**Description :**

Medical High Pressure Vessel

**Technical Specification :**

The high-pressure container has good insulation, aluminum alloy is durable, and liquid nitrogen can hold for a long time. The container can store long biological specimens, preserve animals, semen, and biomedical vaccines.

With 6 filter tanks for biological sample storage, it is only suitable for storing liquid nitrogen and not suitable for transportation.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265