

Product Name :
Medical Blood Bank Refrigerator

Product Code :
ASTRALAB-SUPPLYA13003



Description :

Medical Blood Bank Refrigerator

Technical Specification :

High-brightness digital temperature display, accurate temperature display up to 0.1°C.

The lockable door can prevent others from opening the door without authorization.

It has aluminum thin copper pipes and efficient air-cooled evaporator.

All-round alarm system, with high/low temperature alarm, power failure alarm, access control alarm, power failure alarm.

The two-layer glass window door with automatic defrost function ensures uniform temperature.

Non-CFC rigid polyurethane insulation materials can avoid excessive heat generation.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265