

Product Name :
Medical Vertical Autoclave

Product Code :
ASTRALAB-SUPPLYA15003



Description :

Medical Vertical Autoclave

Technical Specification :

The vertical autoclave is composed of a heating system, a microcomputer control system, and an overheating and overpressure protection system. It has the characteristics of reliable sterilization effect, convenient operation and energy saving. They are ideal equipment for disinfection of surgical instruments, fabrics, glasses, culture media, etc. in clinical, scientific research institutions and other tissues.

All stainless steel structure.

Automatically discharge cold air, and automatically discharge steam after sterilization.

After disinfection, it will automatically shut down and beep.

Quick-opening door structure of hand wheel type.

Automatic protection from overheating and overvoltage.

Easy to operate, safe and reliable.

Door safety lock system.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265