

**Product Name :**  
Astra Pipet Fillers

**Product Code :**  
ASTRALAB-SUPPLYA18004



**Description :**

Astra Pipet Fillers

**Technical Specification :**

Compatible with most of the plastic and glass pipettes from 0.1 -100mL

Six speeds of selection for aspiration and dispensing different liquids

Large LCD display showing low battery warning and speed settings

Enables single-handed operation with minimum effort

Light and ergonomic design providing easy usability

High capacity Li-ion battery enable long runtime of operation

Powerful pump fills a 25mL pipette in <5 seconds

0.45µm replaceable hydrophobic filter

Rechargeable during use



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265