

Product Name :
Auto Portable Low Rise Lifts

Product Code :
ASTRALAB-SUPPLYA20006



Description :

Auto Portable Low Rise Lifts

Technical Specification :

5 position mechanical safety lock system

Portable trolley supports pump and moves lift

Easily adjustable sliding radius arms

Engagement release mode: Automatic mechanical release by novel trip-over mechanism at the highest position

Maximum Load(tons): 4-6

Application: Servicing

Maximum Rise(m): 0-5

After Sales Service: Provided

Max Lifting Weight: 4 ton

Power Source: single phase

Lifting Height: 4.75 feet

Power Of Motor: 2kw

Body Material: steel

Capacity: 4000kg

Min Height: 118 mm

Max Lifting Height: 1152 mm

Max Lifting Height With Adapters: 1260 mm

Overall Length: 2230 mm

Overall Width: 1018 mm

Lifting Time: 35 Sec

Lowering Time: 25 Sec



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265