

**Product Name :**  
Auto Hydraulic Body Jack

**Product Code :**  
ASTRALAB-SUPPLYA20017



**Description :**

Auto Hydraulic Body Jack

**Technical Specification :**

Capacity: 5 to 20 ton, 11-40 Ton

Weight (kg): 30kh

Power (W): no

Color: red

Material: Stainless Steel

Working Pressure: Max. 700 Bar

Hydraulic body jack is not lifting jack but it work with push and pull function which is achieved with a hydraulic pump . it comes with an configuration of 10 ton and 20 ton capacity, all this feature make this tool to have in any body shop.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265