

**Product Name :**  
Auto Mini Portable Tilt Lift

**Product Code :**  
ASTRALAB-SUPPLYA20020



**Description :**

Auto Mini Portable Tilt Lift

**Technical Specification :**

Easy to assemble and ready to use lift.

Revolutionary invention in lifting technology.

Low on budget

No concrete foundation needed like two post lift

Easy to transport, low cost on delivery

Labor saving

Lift can be moved any where in the workshop and a bay can be created.

Maintenance free since no complicated circuits switches etc.

Great saving and value for money.

Can be lifted manually or with help of a drill machine.

Maximum Load(tons): 2-4

Operating Height: 10 feet

Maximum Rise(m): 5-10

Body Material: Stainless Steel

Application: Servicing



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265