

Product Name :
Auto Hydraulic Jack

Product Code :
ASTRALAB-SUPPLYA20023



Description :

Auto Hydraulic Jack

Technical Specification :

Heavy gauge steel side plates prevent twisting

Front wheels with roller bearings and ball bearing mounted rear castor wheels provide stability and easy maneuverability.

Heavy duty pumps, rams, and cylinders are precision machined and inspected to ensure long service life and superior performance

Overload and bypass safety valves prevent the jack from being used beyond its rated capacity and ram over extension

Foot pump pedal raises the lift arm to the load quickly and "T" handle pumps jack under load

(tons) Low Height

(mm) Lift

(mm) Max.Height

(mm) N.W

(KG) Package size

(cm) QTY/CTN

(PCS) G.W

(KG)

2	140	660	800	62	126x44x21	1	66
3	145	465	630	63	126x39x23	1	67
5	160	400	560	97	156x40x28.5	1	107
10	170	400	570	140	169x50x335	1	151
20	190	430	20	295	175x52x35	1	

Capacity 2-3 ton

Low Height 145 mm

Lift 660 mm

Max Height 800 mm

New Weight 62 kg

Package Size 126x44x21 cm

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

