

**Product Name :**  
Auto Steam Vacuum Cleaner

**Product Code :**  
ASTRALAB-SUPPLYA20031



**Description :**

Auto Steam Vacuum Cleaner

**Technical Specification :**

Air Aromatization: Yes

Voltage (V): 230

Weight (kg): 20kg

Automation Grade: Automatic

Boiler Capacity (L): 2

Pressure (Bar): 5

Boiler Pressure (bar): 5

Tank Capacity (l): 2

Boiler Temperature (deg. Celsius): 190

Frequency (Hz): 50

Temperature (deg. Celsius): 2150

Cleaning Type: Steam

Boiler Power (W): 3kw

Number Of Wheels: 4

Cable Length (m): 5

Power Source: Electric

Color: Silver

Temperature: 148 Degree C

Power: 1400W

steam Vacuum cleaner is a professional steam generator and generates a pressure of 5 bar of steam in continuous fashion. this steam cleaner generates steam at 150deg centigrade. the machine also has a feature of Vacuum with heavy duty Vacuum motor which can be used with water curtain as a filter where all the dust can be collected. so you can actually steam and Vacuum at the same time with the help of different adapters provided with machine.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265