

**Product Name :**  
Sintering Furnace

**Product Code :**  
ASTRALAB-SUPPLYA21065



**Description :**

Sintering Furnace

**Technical Specification :**

Some of its advantages include supporting superior thermal treatment of powder, supporting needs of enhancing strength through achieving better bonding together of particles, optimum control support for options related to heating rate, temperature setting and other parameters to achieve reproducible results and others. translucency, thermal conductivity and electrical conductivity properties. The precision process support provided by the system allows the reduction in issues related to porosity as well as enhancing its overall strength,

Sintering Furnace offered can work with maximum temperatures of 1300oC and perfectly supports sintering process involving the formation of solid material mass through heat without melting the material to point of liquification.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265