


**Product Name :**  
Creep Testing Machine

**Product Code :**  
ASTRALAB-SUPPLYA21084



**Description :**

Creep Testing Machine

**Technical Specification :**

Flat specimens made of lead and PE. Tensile stress range 5...25N/mm<sup>2</sup>. Table unit for the investigation of the creep behaviour of different materials. Specimen cross-section 2x5mm, measured length 25mm. Specimen cross-section 2x5mm, measured length 25mm. Position measurement gauge 0...10mm, resolution 0.01mm.

Experiments above or below room temperature possible with conditioning box

Tensile stress range: 5...25N/mm<sup>2</sup>.

Position Measurement.

Measured length: 25mm.

Resolution 0.01mm.

Range: 0...10mm.

Material: lead, PE.

Cross-section: 2x5mm.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265